



The Vicair® Academy AllRounder offers ultimate freedom to the active wheelchair user!

*the best
in
body
support*

**Academy
AllRounder**

Vicair® Technology



The Vicair® technology uses air-filled SmartCells™, moving easily against each other. This way they conform precisely to the body's contours. The individual deformation of the SmartCells™ results in stable body support, combined with superior pressure distribution. The AllRounder consists of a five compartment Academy cushion which can be buckled on and which supports and protects the buttocks and thighs against the effects of sitting on a hard surface, and against impact and abrasion. An upward extension at the rear, containing a layer of tough impact-absorbing foam, provides additional support and protection to the lower back.

Sizes

The AllRounder comes in just one standard size.

Outer cover

This cover, specially designed for the AllRounder, is liquid proof and stretches in two directions.

Distributed by:

VICAIR®

Vicair BV
Koetsierstraat 14
1531 NX Wormer
The Netherlands
T +31(0)75 642 9992
F +31(0)75 642 9998
E info@vicair.com
I www.vicair.com

Usage

Experience unlimited freedom - safely! Don't let the lack of skin protection put you off doing the things you want to do! There's now a product you can rely on to give you the freedom to do all those things you always wanted to. The Vicair® Academy AllRounder has been developed for



you to use whenever you are out of your wheelchair and as a replacement cushion for outdoor activities. The AllRounder enables you to move around the house without your wheelchair. For example, you can get down on the floor to play with your children or to fix the car. The AllRounder protects you not only whilst driving or during a visit to the dentist, but also while you are sailing, canoeing, horseback riding or taking part in any other sporting activities. Don't hold back - do whatever you feel like doing!

Maintenance

The AllRounder's inner cover can easily be cleaned with a standard disinfectant or detergent. The outer cover can be machine washed at 60°C.

Warranty

The Vicair® Academy AllRounder comes with a two year warranty against defects in workmanship and materials. See user manual for more details.

Vicair's other products

