

wheelchairs for function and independence

# Astrotilt<sup>\*</sup>

by PDG

## *tilt in space*

*Wheelchair*

a lightweight, functional, folding  
adult tilt chair



Astrotilt, adult folding chair



50° tilt  
-5° to 45°



Access to rear  
wheels in tilt



Folding or Rigid  
Frame

\*Patent pending.

# Astrotilt<sup>\*</sup>

by PDG

## Features

- **Tilt-In-Space** - 50° for positioning and pressure relief
- **Growable** - maintain a great fit for years
- **Enhances Independence** - for excellent mobility

### Multi-adjustable

#### *Wheelchair*

- Adjustable seat height, seat width, seat depth, back angle, and rear wheel

### Maneuverable

#### *Access*

- Short wheel base

### Lightweight

#### *Tilt*

- Fits in the trunk

### Growable

#### *Frame*

- Ensures long-term fit

### Transport Ready

#### *Brackets*

- Approved tie-down system

## tgif<sup>Warranty</sup> —

- **PDG's exclusive TGIF—Thank Goodness it Fits—Warranty** – ensures our product fits the user
- **Frame Guarantee** – frame has a five year guarantee against defects in materials and workmanship

Seat Width (adjustable)	12" - 18"
Seat Depth (adjustable)	12" - 18"
Seat Height (adjustable)	12" - 20"
Back Angle (adjustable)	84° - 114°
Back Height	14" - 22"
Dynamic Tilt	-5° - +45°
Overall Width	Seat Width + 10"
Overall length (without front rigging)	22"
Front Rigging	70°, 80°, 90° and knee angle adjustable
Transit Options	Available
Weight Capacity	200 lbs (including seating)
Chair Weight <sup>1</sup>	32 lbs

<sup>1</sup>without rear wheels, front rigging, armrests, anti-tippers  
Custom Options Available. See Current Specification Sheets at [www.PDGmobility.com](http://www.PDGmobility.com) for details.



### Product Design Group Inc.

TOLL FREE **1.888.858.4422**  
TOLL FREE FAX **1.877.323.9093**  
[www.PDGmobility.com](http://www.PDGmobility.com) [info@PDGmobility.com](mailto:info@PDGmobility.com)

USA PMB #173, 1124 Fir Avenue Blaine, WA USA 98230  
CANADA 102 - 366 East Kent Ave South Vancouver, BC Canada V5X 4N6

PDG recommends all wheelchairs be fitted to the consumer under the guidance of a health care professional, knowledgeable of the particular needs of the consumer and the limitations of the product.