

wheelchairs for function and independence

Astrotilt Jr.^{*}

by PDG

tilt in space Wheelchair —

brings together all the essential ingredients needed to provide support and versatility for children with changing needs.



50° of tilt
-5° to 45°



Access a table
or desk
while tilted



Folding Frame,
converts to a
rigid chair



Astrotilt Jr., folding

Astrotilt Jr.*

by PDG

Features

Multi-adjustable Wheelchair

- Adjustable seat height, seat width, seat depth, back angle, and rear wheel position

Maneuverable

- Easy maneuverability
- 15" minimum wheel base makes all the difference

Attendant Wheel Locks

- Foot operated

Growable Frame

- Ensures long-term fit

Transport Ready Brackets

- ANSI / RESNA / ISO approved tie down system

Quick-recline Back

- Offers dual position for alternate activities

- **Tilt-In-Space** - 50° for positioning and pressure relief
- **Grows With Your Child** - maintain a great fit for years
- **Fits In Small Spaces** - for easy transport
- **Enhances Independence** - for ease of mobility

tgif Warranty —

- **PDG's exclusive TGIF—Thank Goodness it Fits—Warranty** – ensures our product fits the user
- **Frame Guarantee** – frame has a five year guarantee against defects in materials and workmanship

Seat Width (adjustable)	12" - 18"
Seat Depth (adjustable)	12" - 18"
Seat Height (adjustable)	12" - 20"
Back Angle (adjustable)	84° - 114°
Back Height	14" - 22"
Dynamic Tilt	-5° - +45°
Overall Width	
12" and 16" wheels	Seat Width + 9"
20" and 22" wheels	Seat Width + 10.5"
Overall length (without front rigging)	22"
Front Rigging	70°, 80°, 90°
Transit Options	Available
Weight Capacity	200 lbs (including seating)
Chair Weight†	32 lbs

†without rear wheels, front rigging, armrests, anti-tippers
Custom Options Available. See Current Specification Sheets at www.PDGmobility.com for details.



Product Design Group Inc.

TOLL FREE **1.888.858.4422**
TOLL FREE FAX **1.877.323.9093**
www.PDGmobility.com info@PDGmobility.com

USA PMB #173, 1124 Fir Avenue Blaine, WA USA 98230
CANADA 102 - 366 East Kent Ave South Vancouver, BC Canada V5X 4N6

PDG recommends all wheelchairs be fitted to the consumer under the guidance of a health care professional, knowledgeable of the particular needs of the consumer and the limitations of the product.