

wheelchairs for function and independence

# Eclipse 600

by PDG

## extra-wide, bariatric

*Wheelchair For —*

heavier than average individuals. The Eclipse carries extra heavy loads while delivering style and performance.



Forward wheelbase



Transportable



Narrow Overall Width





## 1000 Pound

*Capacity*

- For extra wide applications

# Eclipse 600

by PDG

## Features

- **Unique Wheelbase Design** – ensures most of the weight is carried on rear wheels, allowing front casters to pivot and roll with ease
- **Rigid Frame** – provides high strength stiffness and long life
- **Easily Transportable** – with removable back, rear wheels, armrests and front rigging
- **Narrow Overall Width** – facilitates accessibility in the environment

## Multi-adjustable

*Wheel chair*

- Features adjustable seat height, seat depth, seat width, back angle and rear wheel

*tgif* Warranty —

- **PDG's exclusive TGIF—Thank Goodness it Fits—Warranty** – ensures our product fits the user
- **Frame Guarantee** – frame has a five year guarantee against defects in materials and workmanship

## Reclining back

*Option*

- Optional reclining back, can be retrofit

	Eclipse 350	Eclipse 600
Seat Width (adjustable)	16" - 24"	22" - 48"
Seat Depth (adjustable)	16" - 20"	16" - 24"
Seat Height (adjustable)	14" - 20"	14" - 20"
Back Angle (adjustable)	85° - 100°	85° - 100°
Recliner	Reclines to 130°	Reclines to 130°
Back Height	18" - 25"	18" - 25"
Overall width	Seat width + 7"	Seat width + 7"
Overall length	25.5" - 28.5"	31.5" - 34.5"
Front Rigging	70° Swing-away, elevating, center mount	60° Swing-away, elevating, center mount
Weight Capacity	350 lbs	1000 lbs
Chair Weight†	32 lbs	35 lbs

†without backposts, armrests, front rigging, rear wheels, anti-tippers

## Width

*Adjustable*

- For long-term fit

## Front-rigging

*Options*

- One-piece, center-mount, swing-away, elevating, for optimal foot positioning

## Solar cushions

- Specially designed for extra wide applications



### Product Design Group Inc.

TOLL FREE **1.888.858.4422**

TOLL FREE FAX **1.877.323.9093**

[www.PDGmobility.com](http://www.PDGmobility.com) [info@PDGmobility.com](mailto:info@PDGmobility.com)

USA

PMB #173, 1124 Fir Avenue Blaine, WA USA 98230

CANADA

102 - 366 East Kent Ave South Vancouver, BC Canada V5X 4N6