

wheelchairs for function and independence

Eclipse 350

by PDG

extra-wide

Wheelchair For —

heavier than average individuals.

The Eclipse 350 carries extra heavy loads while delivering style and performance.



Forward wheelbase



Transportable



Narrow Overall Width



Eclipse 350 by PDG Features



600 Pound
Capacity

- With spoke wheel



350 Pound
Capacity

- With composite wheel



Adjustable
Wheel chair

- Features adjustable seat depth, back angle



Widths
Available

- 18", 20", 22", 24", 26"



Solar cushions

- Specially designed for extra wide applications

- **Unique Wheelbase Design** – ensures most of the weight is carried on rear wheels, allowing front casters to pivot and roll with ease
- **Rigid Frame** – provides high strength stiffness and long life
- **Easily Transportable** – with removable back, rear wheels, armrests and front rigging
- **Narrow Overall Width** – facilitates accessibility in the environment

tgif Warranty —

- **PDG's exclusive TGIF—Thank Goodness it Fits—Warranty** – ensures our product fits the user
- **Frame Guarantee** – frame has a five year guarantee against defects in materials and workmanship

Seat Width	18", 20", 22", 24", 26"
Seat Depth (adjustable)	16" - 20"
Seat Height	14" or 17"
Back Angle (adjustable)	85° - 100°
Back Height	19" or 24"
Overall width	Seat width + 7"
Overall length	31"
Front Rigging	Plug-in, 70° Swing-away or elevating
Weight Capacity	
Composite Wheel	350 lbs
Spoke Wheel	600 lbs
Chair Weight*	34 lbs

*without backposts, armrests, front rigging, rear wheels, anti-tippers



Product Design Group Inc.

TOLL FREE **1.888.858.4422**

TOLL FREE FAX **1.877.323.9093**

www.PDGmobility.com info@PDGmobility.com

USA PMB #173, 1124 Fir Avenue Blaine, WA USA 98230

CANADA 102 - 366 East Kent Ave South Vancouver, BC Canada V5X 4N6

PDG recommends all wheelchairs be fitted to the consumer under the guidance of a health care professional, knowledgeable of the particular needs of the consumer and the limitations of the product.