

# MaTRx seating series:

## matrx-v

## Kid\*ab\*ra



Vahen, an Administrative Coordinator and part-time model, was diagnosed with Transverse Myelitis in 1999.



The popular MaTRx-V and Kid\*ab\*ra cushions and bases are carefully designed using a combination of HR and highly sensitive visco-elastic foam. They provide extremely effective pressure reduction and positioning — without any of the hassles of maintenance or adjustment.

## MaTRx-V/Kid\*ab\*ra cushion/base

*standard features & options*

### **VISCO LAYER**

Highly sensitive molded visco-elastic foam top layer

- ▶ Provides exceptional pressure reduction by maximizing overall surface contact

### **BOTTOM LAYER**

Anatomically contoured, lightweight, molded HR polyurethane foam base

- ▶ Anatomically contoured to promote good positioning and stability

### **INNER COVER**

Durable inner cover protects foam

- ▶ Zippered to allow cushion modifications

### **OUTER COVER**

Specially designed three-layer outer cover

- ▶ Moisture proof/breathable to keep clients cool and dry
- ▶ Kid\*ab\*ra cover is custom printed with appealing design that kids will love

### **BASE (OPTIONAL)**

- ▶ MaTRx-V/Kid\*ab\*ra cushions are available with a built-in drop base and choice of drop hooks

### **PELVIC OBLIQUITY KIT (POK)**

- ▶ Easily installed in all MaTRx-V and Kid\*ab\*ra cushions/bases
- ▶ Accommodates 1/2" or 1" obliquities

### **SEAT DEPTH REDUCTION/LEG LENGTH DISCREPANCY**

- ▶ Foam, inner and outer covers modified to your specifications upon request

### **SIZES**

- ▶ Standard sizes from 12" - 22" wide, 12" - 20" deep. See price list for more information.



	Kid*ab*ra Cushion	Kid*ab*ra Base	MaTRx-V Cushion	MaTRx-V Base
Weight capacity	150 lb	150 lb	250 lb/300 lb	250 lb
Cushion width	12", 14"	12", 14"	14", 16", 18", 20", 22"	16", 17", 18", 19", 20"
Cushion depth	12", 14"	12", 14"	15", 16", 17", 18", 20"	16", 18"
Cushion height (at leg trough)	1½"	1½"	2"	2"
Includes (built-in drop base & hooks)	—	●	—	●

## Options

