

MaTRx seating series:

matrx-pb



NEW!



MaTRx-PB Posture back & Elite back *standard features & options*

MaTRx PB POSTURE BACK

- ▶ Comfortable contour provides excellent midline support and positioning
- ▶ Available in Standard (5") or Deep (7") contour depth
- ▶ Standard sizes 14" - 20" wide, 12" - 20" high
- ▶ Heavy Duty sizes 20" - 24" wide, 16" - 20" high
- ▶ 300/500 lb weight capacity

MaTRx PB ELITE BACK

- ▶ Sleek, stylish and breathable back shell
- ▶ 3" contour depth for more active lifestyles
- ▶ Optional formable, width adjustable lateral supports
- ▶ Standard sizes 12"-20" wide, 12"-20" high
- ▶ 300 lb weight capacity

"EASY SET HARDWARE"

Standard on all MaTRx Posture Backs/Elite Backs
Fast & easy to install and adjust

- ▶ "Infinite" adjustment of Recline (up to 35°) and seat depth (up to 4") and even accommodates rotation!
- ▶ Width adjustable mounting pins allow adjustment to fit up to 3 chair widths
- ▶ Quick release hardware facilitates portability
- ▶ All adjustments can be made with client in chair

MaTRx Contour back

- ▶ Comfortable and simple
- ▶ Mildly contoured posterior back support
- ▶ Fixed hardware ensures back remains securely on wheelchair
- ▶ Limited depth/Recline adjustability
- ▶ Sizes from 15" - 22" wide, 12" - 17" high

The MaTRx backs take postural control and comfort to a whole new level. Designed for everyone from active users to more involved clients, from pediatric through geriatric, we have a back model to meet the widest range of needs.

The MaTRx Contour Back offers comfortable simplicity with fixed mounting hardware, while the MaTRx PB Posture Backs and Elite Backs feature postural support and comfort, along with our highly popular "Easy Set" hardware for infinite adjustment, portability and ease of installation.

MaTRx PB Posture Back/Elite Back hardware

All adjustments can be made with user seated in w/c — only one tool required!



"Easy Set" hardware can be installed and adjusted in minutes

'Infinite' adjustment of back angle (35°), depth (4"), and even rotation!

Unique, adjustable mounting pins allow each back to fit chairs one inch wider or narrower!

2 year warranty

Hardware only available on MaTRx Series backs



Kid*ab*ra fabric optional on 14" widths



	MaTRx-PB & PB Deep Backs	MaTRx-PB Heavy Duty & Heavy Duty Deep Backs	MaTRx-PB Elite Back/Elite with Laterals	Contour Back
Weight capacity	300 lb	500 lb	300 lb	300 lb
Back width	14", 16", 18", 20"	20", 22", 24"	12" - 20"	15" - 22"
Fits w/c widths	13" to 20"	19" to 25"	12" - 20"	15" - 22"
Back height	12", 16", 20"	16", 20"	12", 16", 20"	12" - 17"
Contour depth	5" / 7"	5" / 7"	3"	2"

Options

Available on MaTRx PB Posture Backs, Elite Backs and Contour Backs

Extended Shear Reduction Option (ESR)

Easy to use kit allows customization of back shape. Available on MaTRx PB Posture Backs only.

Foam in Place Option (FIP)